

Ratchets

Refer to p: 38-58

MIS offers a line of uniquely designed ratchets, to simplify both prosthetic screw tightening and implant insertion, allowing an accurate and safe performance. MIS ratchets are high durability instruments with a unique mechanism that simplifies use and cleaning.

Torque Ratchet for prosthetic screws

MT-RI040



Controlled tightening torque of prosthetic screws (10-30Ncm).

Torque Ratchet for implant placement

MT-RT070



Controlled tightening torque of implant placement (30-75 Ncm). The recommended torque for implant placement is 35-60 Ncm.

Monoblock Ratchet

MT-RI030



The monoblock ratchet can be used for implant placement and tightening or loosening screws.